

Identification Reference Chart



1 Ask for part number - Part numbers are located on all product labels, patient guides and brochures

<p>Aspira* Drainage System</p>	<p>Aspira* Drainage System with Blue Valve</p>
---------------------------------------	---

2 Product Code/Reorder Number

<p>Aspira* Drainage System</p>	<p>Aspira* Drainage System with Blue Valve</p>
---------------------------------------	---

3 Drainage Cap

<p>Aspira* Drainage System</p>	<p>Aspira* Drainage System with Blue Valve</p>
---------------------------------------	---

NOTE: Steps 3 & 4: When the Aspira* Drainage System's first generation catheter is repaired with the Aspira* Drainage System's second generation "Blue Valve" repair kit (Product Code 4992306), the Aspira* Drainage System's second generation "Blue Valve" compatible accessories must be ordered.

4 Ask customer/patient to remove the dressing kit bandage and ask them what color the valve is?

<p>Aspira* Drainage System</p>	<p>Aspira* Drainage System with Blue Valve</p>
---------------------------------------	---

5 Ask customer/patient to remove the dressing kit bandage and ask them if there is printing on the catheter.

<p>Aspira* Drainage System</p>	<p>Aspira* Drainage System with Blue Valve</p>
---------------------------------------	---

4991502 - (5 Qty) 1,000 mL Aspira* Drainage Bag	4992301 - (5 Qty) 1,000 mL Aspira* Drainage Bag
4991503 - Dressing Kit - 5 Kits included	4991503 - Dressing Kit - 5 Kits included

NOTE: Both systems use the same dressing kit



Bard Access Systems, Inc.
 605 North 5600 West
 Salt Lake City, UT 84116
 801-522-5000, 800-545-0890
 www.bardaccess.com

* Bard and Aspira are trademarks and/or registered trademarks of C. R. Bard, Inc. © 2011 C. R. Bard, Inc. All rights reserved.